

COMPOSITE MICROFILL POLISHERS - SINGLE USE


# DISPOSABLE KENDA COMPOSITE MICROFILL 2-step polishing system

For single patient use.




Top selling polisher for all microfilled composite restorations. Achieves an optimal and long lasting high shine.

Also, suited for simply removing cement remnants of orthodontic bands (brackets).

Disposable polisher saves sterilisation cost and time.



Order REF	PS 4009	PS 4011	PS 4012	PS 5009	PS 5011	PS 5012
Size Ø x L mm	10x1	6x7,5	5x10	10x1	6x7,5	5x10
ISO No.	638 204... 303 512 100	030 512 060	243 512 050	303 492 100	030 492 060	243 492 050
Grit	fine			ultra fine		
Color	yellow			blue		
⌚	5000 RPM			5000 RPM		

-  Box of: 6, 12, 25, 100
-  Pressure: 0,3 - 1 N
- Shank: RA, medical plastic
- Medical device: class IIa
-  disposable