

POLISHERS FOR FLEXIBLE DENTURES


# KENDA RED

## 1-step polishing system

For adjusting and polishing dentures and flexible partials on chrome cobalt alloys without damaging the alloy, mounted and unmounted.

Recommendation:

For high shine use KENDA 244 BLUE polishing paste and slimline brush (page K-39)



Order REF	RED 8	RED 14	10008	10014	301-L
Size Ø x L mm	12,5x22	10x24.5	12,5x22	10x24,5	5x13
ISO No.	658 104 ... 251 525 125	497 525 100	658 900 ... 251 525 125	497 525 100	312 104 ... 610 415 050
Grit	medium		medium		
Color	red		red		
⊙	10'000 - 15'000 RPM		10'000 RPM		



Box of: 6, 12, 25, 100  
Pressure: 0,3 - 2 N  
Shank: HP, stainless steel