

KENDA ZIRCO1

Only **1** step to a high luster finish on zirconia and all ceramics

new



- ✓ *First of its kind: Finally achieve exceptional results on zirconia, lithium disilicate and feldspathic porcelain in just one fast step.*
- ✓ *Use variable pressure techniques, familiar from 1-step composite polishing systems, to transition from scratched surfaces to a high-shine finish.*
 - More pressure @ 15'000rpm: Smoothing
 - Light touch @ 7'500rpm: High Shine Polishing
- ✓ *Produces very low surface roughness values and highly esthetic gloss values.*
- ✓ *Decreases plaque accumulation, reduces chipping and fracture of restorations and prevents excessive wear of opposing enamel.*
- ✓ *Available in the two most common shapes for polishing intraoral adjustments on ceramic restorations.*



REF 0233



REF 0240



Set REF 0230

KENDA AG

Dental Manufacturing · Industriezone Neugut
Kanalstrasse 18 · LI-9490 Vaduz · Principality of Liechtenstein

Phone +423 388 23 11 · Fax +423 388 23 19
info@kenda-dental.com · www.kenda-dental.com

CE
0120

Edition 03/2017